

Sorenson OneTouch™ Tips

Bring greater comfort to your pipetting experience and reduce wear on your equipment with versatile OneTouch pipet tips.

- ▶ **Ergonomic:** Soft top collar reduces the sealing and ejection force required, making pipetting more comfortable and efficient
- ▶ **Accurate:** Rigid, ultra-clear polypropylene body provides exceptional clarity for sample visualization and accurate measurements
- ▶ **Versatile:** Wide compatibility means you can employ them with various pipettor models, allowing for seamless transitions between equipment
- ▶ **Easy identification:** Pipet tips are color-coded by volume, allowing for quick and error-free identification
- ▶ **Long-lasting pipettor performance:** OneTouch tips are engineered to minimize wear on pipettor instruments, ultimately extending their lifespan and reducing maintenance needs
- ▶ **Certified quality:** The tips are certified RNase-/DNase-free
- ▶ **Variety of options:** Choose from barrier or standard tips in various sizes to suit your specific needs

Ordering Information

Sorenson OneTouch Tips

● 10 µL OneTouch Tips

Cat. No.	Description	Packaging	Tips/Pack	Packs/Case
10300	10 µL OneTouch tip	Cell rack	960	5
10310	10 µL OneTouch tip, sterile	Cell rack	960	5
10320	10 µL OneTouch tip, sterile, barrier	Cell rack	960	5

● 200 µL OneTouch Tips

Cat. No.	Description	Packaging	Tips/Pack	Packs/Case
10330	200 µL OneTouch tip	Cell rack	960	5
10340	200 µL OneTouch tip, sterile	Cell rack	960	5
10350	200 µL OneTouch tip, sterile, barrier	Cell rack	960	5

● 1000 µL OneTouch Tips

Cat. No.	Description	Packaging	Tips/Pack	Packs/Case
10360	1000 µL OneTouch tip	Hinged rack	768	5
10370	1000 µL OneTouch tip, sterile	Hinged rack	768	5
10380	1000 µL OneTouch tip, sterile, barrier	Hinged rack	768	5

