

Sorenson[®] Pipette Tips

Innovative, ergonomic, environmentally friendly





Sorenson® Pipette Tips

Sorenson MμLtiRack Tip System	2
Racked, Filtered, Low Binding	2
Racked	2
TipStation	3
Reload Stack	3
Sorenson HybridRack™ Pipette Tips	4
Racked, Filtered, Low Binding	5
Racked	5
Sorenson LTS Compatible Pipette Tips	6
Racked, Filtered, Sterile	7
Racked, Sterile	7
Racked	7
Sorenson HybridRack, Filtered, Sterile	7
TipStation, Sterile	7
TipStation	7
Reload Stack	7

Sorenson® MμLtiRack Tip System

Engineered to fit most pipettors on the market, the Sorenson MμLtiRack Tip System is a complete pipette tip system offering in various packaging styles to meet every lab's needs. The MμLtiRack Tip System offers standard and extended-length sizes in volumes from 0.1 μL to 1250 μL.

Racked and reload packaging styles feature preloaded tips in easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and maximizes lab space. Sorenson MμLtiRack Tip System products are made in an ISO 13485 and ISO 14001 registered facility.

- ▶ Fine point
- ▶ Graduated
- ▶ RNase-/DNase-free
- ▶ Human DNA-free
- ▶ PCR inhibitor-free
- ▶ Non-pyrogenic

Sorenson MμLtiRack Tip System

Racked, Filtered, Low Binding



Cat. No.	Volume	Extended Length	Sterile	Low Binding	Quantity
77250T	10 μL	Yes	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77260T	20 μL	No	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77270T	100 μL	No	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77280T	200 μL	No	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77290T	200 μL	Yes	Yes	Yes	96/rack, 8 racks/pack, 5 packs/case
77310T	1250 μL	Yes	Yes	Yes	96/rack, 8 racks/pack, 5 packs/case

Racked



Cat. No.	Volume	Extended Length	Sterile	Low Binding	Quantity
77840	10 μL	Yes	No	No	96/rack, 10 racks/pack, 5 packs/case
77850	10 μL	Yes	Yes	No	96/rack, 10 racks/pack, 5 packs/case
77850T	10 μL	Yes	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77860	200 μL	No	No	No	96/rack, 10 racks/pack, 5 packs/case
77870	200 μL	No	Yes	No	96/rack, 10 racks/pack, 5 packs/case
77870T	200 μL	Yes	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77880	300 μL	No	No	No	96/rack, 10 racks/pack, 5 packs/case
77890	300 μL	No	Yes	No	96/rack, 10 racks/pack, 5 packs/case
77890T	300 μL	No	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77300T	300 μL	Yes	Yes	Yes	96/rack, 8 racks/pack, 5 packs/case
77900	1250 μL	Yes	No	No	96/rack, 8 racks/pack, 5 packs/case
77910	1250 μL	Yes	Yes	No	96/rack, 8 racks/pack, 5 packs/case

TipStation



Cat. No.	Volume	Extended Length	Sterile	Low Binding	Quantity
27780	10 μL	Yes	No	No	192/layer, 5 layers/pack, 5 packs/case
27790	10 μL	Yes	Yes	No	192/layer, 5 layers/pack, 5 packs/case
27760	200 μL	No	No	No	192/layer, 5 layers/pack, 5 packs/case
27770	200 μL	No	Yes	No	192/layer, 5 layers/pack, 5 packs/case

Reload Stack



Cat. No.	Volume	Extended Length	Sterile	Low Binding	Quantity
70500	10 μL	Yes	No	No	96/deck, 10 decks/pack, 4 packs/case
70510	200 μL	No	No	No	96/deck, 10 decks/pack, 4 packs/case
70530	300 μL	No	No	No	96/deck, 10 decks/pack, 4 packs/case
70540	1250 μL	Yes	No	No	96/deck, 8 decks/pack, 4 packs/case

Sorenson® HybridRack™ Pipette Tips

Have you ever wanted your lab to use less plastic? Have you been looking for an easy way to be more sustainable around the lab?

Just like you, Sorenson takes sustainability seriously and is always looking for ways to reduce our carbon footprint. The HybridRack does just that, marrying functionality and sustainability by pairing a polypropylene deck with a paperboard base to allow for more sustainability in science. All without sacrificing quality, comfort, or reliability.

LESS PLASTIC

Innovative paperboard and plastic hybrid design uses **70% less plastic** than a traditional tip rack.



Paperboard base can be easily recycled using most curbside recycling programs which can significantly reduce the amount of plastic being thrown away in the lab.



The unique hybrid design creates a stable platform allowing end users to pipet directly from the paperboard base – even with a multi-channel instrument.


- Same consistent and high quality performance you have come to expect from Sorenson pipette tips
- Unique design allows for use as a filter tip reload system
- Plastic deck is compatible with Sorenson M μ LtiRack Tip System racks
- Made in an ISO 13485 and ISO 14001 registered facility
- Sterile
- Nonpyrogenic
- RNase-/DNase-free
- Human gDNA- and PCR-inhibition-free
- Available with polyethylene filter
- Available with Sorenson Low Binding Polymer surface to improve accuracy and ensure maximum sample recovery

ACT. The Environmental Impact Factor Label		US
Sorenson™ 10uL HybridRack™ Tip Salt Lake City, Utah, United States SKU 72120		
Environmental Impact Scale Decreasing Environmental Impact		
Manufacturing		
Manufacturing Impact Reduction	1.0	
Renewable Energy Use	No	
Responsible Chemical Management	1.0	
Shipping Impact	1.0	
Product Content	5.0	
Packaging Content	1.0	
User Impact		
Energy Consumption (kWh/day)	N/A	
Water Consumption (gallons/day)	N/A	
Product Lifetime	6.0	
End of Life		
Packaging	1.0	
Product	6.6	
Environmental Impact Factor:		22.6
Label Valid Through:		June 2027
		act.mygreenlab.org



my green lab.
ACT Ecolabel

All HybridRack SKUs have been third-party audited to the My Green Lab® ACT® Ecolabel Standard, the premier eco-nutrition label for laboratory products.





Sorenson® HybridRack™ Pipette Tips

Racked, Filtered, Low Binding

Cat. No.	Description	Quantity	Similar Plastic Rack Cat. No.
87240T	10 µL HybridRack extended length filtered tip, racked, low binding, sterile	960/pk, 5 pks/cs	38000T
87250T	20 µL HybridRack filtered tip, racked, low binding, sterile	960/pk, 5 pks/cs	14210T
87260T	100 µL HybridRack filtered tip, racked, low binding, sterile	960/pk, 5 pks/cs	36060T
87270T	200 µL HybridRack filtered tip, racked, low binding, sterile	960/pk, 5 pks/cs	30550T
87830T	200 µL HybridRack filtered tip, extended length, racked, low binding, sterile	768/pk, 5 pks/cs	17370T-X
87280T	300 µL HybridRack extended length filtered tip, racked, low binding, sterile	768/pk, 5 pks/cs	30510T-X
87290T	1000 µL HybridRack extended length filtered tip, racked, low binding, sterile	768/pk, 5 pks/cs	34000


Racked

Cat. No.	Description	Quantity	Similar Plastic Rack Cat. No.
72120	10 µL HybridRack extended length tip, racked, sterile	960/pk, 5 pks/cs	37660
72130	200 µL HybridRack tip, racked, sterile	960/pk, 5 pks/cs	15670
72140	300 µL HybridRack tip, racked, sterile	960/pk, 5 pks/cs	31790
72150	300 µL HybridRack extended length tip, racked, sterile	768/pk, 5 pks/cs	30930
72160	1000 µL HybridRack extended length tip, racked, sterile	768/pk, 5 pks/cs	34760

Sorenson® LTS Compatible Pipette Tips

The new Sorenson LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch™ pipettors, all with the quality you have come to expect from Sorenson.

- ▶ **Comfort You Can Feel:** LTS Compatible Tips reduce the force required to seal and eject, enabling a streamlined experience from aspiration to dispensing. Avoid slamming with excessive force on the pipettor for a proper seal, eliminate the struggle to eject tips, and minimize hand and wrist strain after extensive pipetting.
- ▶ **Convenience of Choice:** LTS Compatible Tips are available in your favorite packaging formats, including racked, racked/filtered, TipStation, reload stacks, and our new more sustainable solution: Sorenson HybridRack™ Pipette Tips. If you're not using an LTS pipettor, no problem. We provide a complete range of universal pipette tips to choose from.
- ▶ **Quality:** Our United States manufacturing facility is audited by third parties to ensure adherence to regulated manufacturing standards, including ISO 13485 and ISO 14001. Sterile products feature a Sterility Assurance Level (SAL) of 10^{-3} , and all products are made using USP Class VI resin.
- ▶ **Peace of Mind:** Our LTS Compatible pipette tips are RNase-/DNase-free, Human gDNA-free, PCR inhibition-free, and nonpyrogenic so you can pipet stress-free, knowing your samples are safe. Additionally, the tips are autoclavable at 121°C for 15 minutes.



Sorenson® LTS Compatible Pipette Tips



Racked, Filtered, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
87940	20 µL LTS compatible tip, filtered, racked, sterile	960	4,800
87950	200 µL LTS compatible tip, filtered, racked, sterile	960	4,800
87960	1000 µL LTS compatible tip, filtered, racked, sterile	768	3,840

Racked, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
87910	20 µL LTS compatible tip, racked, sterile	960	4,800
87920	200 µL LTS compatible tip, racked, sterile	960	4,800
87930	1000 µL LTS compatible tip, racked, sterile	768	3,840

Racked

Cat. No.	Description	Qty/Pk	Qty/Cs
87880	20 µL LTS compatible tip, racked	960	4,800
87890	200 µL LTS compatible tip, racked	960	4,800
87900	1000 µL LTS compatible tip, racked	768	3,840



Sorenson HybridRack™, Filtered, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
87970	20 µL LTS compatible tip, filtered, Sorenson HybridRack, sterile	960	4,800
87980	200 µL LTS compatible tip, filtered, Sorenson HybridRack, sterile	960	4,800
87990	1000 µL LTS compatible tip, filtered, Sorenson HybridRack, sterile	960	4,800

TipStation, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
88040	20 µL LTS compatible tip, TipStation, sterile	960	4,800
88060	200 µL LTS compatible tip, TipStation, sterile	960	4,800

TipStation

Cat. No.	Description	Qty/Pk	Qty/Cs
88030	20 µL LTS compatible tip, TipStation	960	4,800
88050	200 µL LTS compatible tip, TipStation	960	4,800



Reload Stack

Cat. No.	Description	Qty/Pk	Qty/Cs
88000	20 µL LTS compatible tip, reload stack	960	3,840
88010	200 µL LTS compatible tip, reload stack	960	3,840
88020	1000 µL LTS compatible tip, reload stack	768	3,072

