

## 100-1000 $\mu$ l MultiGuard Pipette Tip



### Attributes-

*Material:* Polypropylene, Natural Color

*Maximum volume:* 1000  $\mu$ l

*Certified RNase/DNase Free*

*Certified PCR Inhibitor Free*

*Certified Human DNA Free*

*Certified Pyrogen Free*


*Length:* 2.848 inches

*Orifice Diameter:* 0.035 inches



### Features & Function-

- Standard 1000 $\mu$ l pipette tip with universal fit on most popular laboratory pipettors.
- Available with low-binding surface technology.
- Highly effective hydrophobic aerosol barrier.
- Various packaging configurations available.



### Compatible Pipettors-

- Gilson Pipetman 1000 $\mu$ l
- Biohit Proline 1000 $\mu$ l
- Eppendorf Research 1000 $\mu$ l
- Eppendorf Reference 1000 $\mu$ l
- Labsystems Finnpipette 1000 $\mu$ l
- Oxford/Nichiryo Benchmate 1000 $\mu$ l

*Universal Fit for Today's Laboratories!*

Cat #	Low Binding Cat #	Description	Packaging
14200	14200T	Rack Sterile, 100-1000 $\mu$ l, Natural	1,000 tips/pack; 1 pack/case

