

1000 XT Multiguard Barrier Tip



Attributes-

Material: Polypropylene, Natural Color

Maximum volume: 1250 µl

Certified RNase/DNase Free

Certified PCR Inhibitor Free

Certified Human DNA Free

Certified Pyrogen Free


Length: 4.0191 inches

Orifice Diameter: 0.034 inches



Features and Function-

- Graduation marks at 200 µl, 500 µl, and 1000 µl.
- Racked on 9 mm centers, multichannel compatible.
- Rack features a hinged lid.
- BioGuard hydrophobic aerosol barrier



Compatible Pipettors-

Gilson Pipetman - Pipetman P1000

BioHit, Proline Mechanical pipettors (adjustable volume) - 100 - 1000 µl

BioHit, mLINE Single Channel Mechanical Pipettors (adjustable volume) - 100 - 1000 µl

Eppendorf, Reference Series 2000 adjustable volume - 100 - 1000 µl

Eppendorf, Research Series 2100 adjustable volume - 100 - 1000 µl

ThermoLabSystems-Finnpipette, Finnpipette Digital Single channel - 200 - 1000µl

ThermoLabSystems-Finnpipette, Finnpipette Digital Single channel - 100 - 1000 µl

Universal Fit for Today's Laboratories!

Cat#	Description	Packaging
34000	50-1250 µl, Standard Sterile	768 tips/pack; 5 packs/case
34000T	50-1250 µl, Low Binding Sterile	768 tips/pack; 5 packs/case

