


Sorenson BioScience Product Specification Sheet

ultraAmp™ Real Time PCR Sealing Film



Attributes

Material: Polyester
Certified RNase/DNase Free
Certified PCR Inhibitor Free
-40° to +120° C temp. range
2 mil thickness



Features & Function

- Ultra clear for Real Time PCR applications
- Thin, easy to peel
- Notched corner prevents peeling off of corner wells
- Easy-tear pull tabs
- Versatile: use for all PCR applications
- Low-static for easy application and handling

Cat #
36590

Description
ultraAmp Real Time PCR Sealing Film

Packaging
100 sheets/box

