

MiniFlex™ Gel Loading Tips



Attributes-

Material: Polypropylene, Natural Color

Maximum volume: 10 μ l

Certified RNase/DNase Free

Certified PCR Inhibitor Free

Certified Human DNA Free

Certified Pyrogen Free

Length: 1.800 inches

Orifice Diameter: 0.015 inches

Length of flat section: 0.375 inches



Features & Function-

- Universal fit for research-grade ultramicro pipettors
- Flat tip ideal for loading vertical PAGE Gels
- 45mm length for access to bottom of narrow vessels
- 2 μ l graduation for quick volume check
- Different outside diameters for different gel thickness
- Load 0.6mm gels, 0.4mm gels, or 0.2mm gels easily

Compatible Pipettors-

Gilson Pipetman 2 μ l, 10 μ l


BioHit Proline 2.5 μ l, 10 μ l

Eppendorf Research 2.5 μ l, 10 μ l

Eppendorf Reference 2.5 μ l, 10 μ l

Labsystems Finnpipette 2 μ l, 10 μ l

Oxford/Nichiryo Benchmate 2 μ l, 10 μ l



The end of the flat versions feature a 0.375 inch long "duckbill" that is either .37mm or .17mm thick!

Universal Fit for Today's Laboratories!

Cat #	Description	Packaging
15100	MiniFlex Round (0.57mm o.d.)	200 tips/pack; 4 packs/case
15110	MiniFlex Round (0.57mm o.d.), sterile	200 tips/pack; 4 packs/case
15120	MiniFlex Flat (0.37mm o.d.)	200 tips/pack; 4 packs/case
15130	MiniFlex Flat (0.37mm o.d.), sterile	200 tips/pack; 4 packs/case
17350	MiniFlex Flat (0.17mm o.d.)	200 tips/pack; 4 packs/case
17360	MiniFlex Flat (0.17mm o.d.), sterile	200 tips/pack; 4 packs/case

