

OneTouch 1000µl Barrier Pipette Tip (50-1250µl)




Attributes

Functional Material: Natural Polypropylene
 Collar Material: TPE
 Maximum volume: 1250µl
 Certified RNase/DNase Free
 Certified PCR Inhibitor Free
 Certified Human DNA Free
 Certified Pyrogen Free
 Length: 4.025 inches
 Orifice Diameter: 0.032 inches
 Barrier: small micron, hydrophobic plug

Features and Function

- Ergonomic dual-plastic construction eases mounting and eject force.
- Extra-long reach for deep vessels.
- Graduation marks at 200µl, 500µl, and 1000µl.
- Available racked on 9mm centers, multichannel compatible.
- Plastic rack features a hinged lid; also available in e-dek, our innovative enviro-friendly, recyclable paperboard rack.



Compatible Pipettors

Gilson Pipetman - Pipetman P1000
 BioHit, Proline Mechanical pipettors (adjustable volume) - 100 - 1000µl
 BioHit, mLINE Single Channel Mechanical Pipettors (adjustable volume) - 100 - 1000µl
 Eppendorf, Reference Series 2000 adjustable volume - 100 - 1000µl
 Eppendorf, Research Series 2100 adjustable volume - 100 - 1000µl
 Universally fits most popular 1000µl pipettors

Cat#	Description	Packaging
------	-------------	-----------

Plastic Racks

10380	Sterile, 1000µl OneTouch Barrier Pipette Tip, natural with blue collar	768 tips/pack; 5 packs/case
-------	--	-----------------------------

