

# Sorenson BioScience Product Specification Sheet

## OneTouch 10 $\mu$ l Barrier Pipette Tip (1-20 $\mu$ l)




### Attributes

Functional Material: Natural Polypropylene  
Collar Material: TPE  
Maximum volume: 20 $\mu$ l  
Certified RNase/DNase Free  
Certified PCR Inhibitor Free  
Certified Human DNA Free  
Certified Pyrogen Free  
Length: 1.785 inches  
Orifice Diameter: 0.015 inches  
Barrier: small micron, hydrophobic plug



### Features and Function

- Ergonomic dual-plastic construction eases mounting and eject force.
- Universal fit- get easy load/eject with any common pipettor.
- Graduation marks at 1 $\mu$ l, 5 $\mu$ l, and 10 $\mu$ l.
- Available in plastic CellRacks or in e-dek, our innovative enviro-friendly, recyclable paperboard rack.



### Compatible Pipettors

Gilson Pipetman 2 $\mu$ l, 10 $\mu$ l  
Biohit Proline 2.5 $\mu$ l, 10 $\mu$ l  
Eppendorf Research 10 $\mu$ l  
Eppendorf Reference 10 $\mu$ l  
Labsystems Finnpipette 2 $\mu$ l, 10 $\mu$ l  
Oxford/Nichiryo Benchmate 2 $\mu$ l, 10 $\mu$ l

Cat#	Description	Packaging
<b>Plastic Racks</b>		
10320	Sterile, 10 $\mu$ l OneTouch Barrier Pipette Tip, natural with red collar	960 tips/pack; 5 packs/case

