

Sorenson BioScience Product Specification Sheet

OneTouch 200 μ l Barrier Pipette Tip (1-330 μ l)




Attributes

Functional Material: Natural Polypropylene
Collar Material: TPE
Maximum volume: 330 μ l
Certified RNase/DNase Free
Certified PCR Inhibitor Free
Certified Human DNA Free
Certified Pyrogen Free
Length: 2.315 inches
Orifice Diameter: 0.018 inches
Barrier: small micron, hydrophobic plug



Features and Function

- Ergonomic dual-plastic construction eases mounting and eject force.
- Universal fit- get easy load/eject with any common pipettor.
- Graduation marks at 50 μ l, 100 μ l, and 200 μ l.
- Available in plastic CellRacks or in e-dek, our innovative enviro-friendly, recyclable paperboard rack.



Compatible Pipettors

Gilson Pipetman 200 μ l
Eppendorf Research 200 μ l
Eppendorf Reference 200 μ l
LabSystems Finnpipette 200 μ l
Oxford/Nichiryo Benchmate 200 μ l

Cat#	Description	Packaging
Plastic Racks		
10350	Sterile, 200 μ l OneTouch Barrier Pipette Tip, natural with yellow collar	960 tips/pack; 5 packs/case

